

General Information

Gene Name Synonym

IgLON family member 4

Protein Construction

A DNA sequence encoding the human NEGR1 (NP_776169.2) (Met 1-Gly 324) without the propeptide was expressed, with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Val 38

Molecular Mass

The secreted mature form of recombinant human NEGR1 consists of 298 amino acids and predicts a molecular mass of 32.9 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rh NEGR1 is approximately 45-48 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human NEGR1 Protein (His Tag) SDS-PAGE