

# H9N2 HA (His Tag) recombinant protein



Catalog Number: 501426

## General Information

### Protein Construction

A DNA sequence encoding the N-terminal segment (Met 1-Arg 338) of the influenza hemagglutinin (Influenza A virus (A/Chicken/Hong Kong/G9/97 (H9N2)) (AAF00701.1), termed as HA1, was fused with a polyhistidine tag at the C-terminus.

### Organism

H9N2

### Expression Host

Human Cells

## QC Testing

### Purity

> 97 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Asp 19

### Molecular Mass

The secreted recombinant HA1 subunit of influenza A H9N2 HA (Influenza A virus (A/Chicken/Hong Kong/G9/97(H9N2))) comprises

331 amino acids and has a predicted molecular mass of 37 kDa. As a result of glycosylation, it migrates as an approximately 50 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Influenza A H9N2 (A/Chicken/Hong Kong/G9/97) Hemagglutinin Protein (HA1 Subunit) (His Tag) SDS-PAGE