

Influenza B HA (B/Massachusetts/03/2010) recombinant protein



Catalog Number: 501431

General Information

Protein Construction

A DNA sequence encoding the full length of Influenza B virus (B/Massachusetts/03/2010) hemagglutinin (AGB07680.1) (Met1-Leu585), termed as HA, was expressed.

Organism

Influenza B

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 16

Molecular Mass

The recombinant hemagglutinin of the Influenza B virus (B/Massachusetts/03/2010) consists 570

amino acids and predicts a molecular mass of 61.7 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 150 mM NaCl, pH 7.5, 5 % glycerol.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza B (B/Massachusetts/03/2010)
Hemagglutinin / HA0 Protein