

# Human Cystatin A / CSTA (His Tag) recombinant protein



Catalog Number: 501510

## General Information

### Gene Name Synonym

Cystatin-AS; Stefin-A; Cystatin-A, N-terminally processed

### Protein Construction

A DNA sequence encoding the human CSTA (P01040) (Ile 2-Phe 98) was expressed with a polyhistidine tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

Please contact us for more information.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Met

### Molecular Mass

The recombinant human CSTA consisting of 108

amino acids and has a calculated molecular mass of 12.4 KDa. It migrates as an approximately 14 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile 50mM Tris, 1mM EDTA, pH 7.5

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human Cystatin A / CSTA Protein (His Tag) SDS-PAGE