

Human CRACC/SLAMF7?CD319? (His Tag) recombinant protein



Catalog Number: 501550

General Information

Gene Name Synonym

CD2 subset 1; CD2-like receptor-activating cytotoxic cells; Membrane protein FOAP-12; Novel Ly9; Protein 19A

Protein Construction

A DNA sequence encoding the human SLAMF7 (NP_067004.3) extracellular domain (Met 1-Met 226) was expressed, fused with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Activity

Measured by its ability to bind biotinylated human SLAMF7 in a functional ELISA.

Purity

> 93 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ser 23

Molecular Mass

The recombinant human SLAMF7 consists of 215 amino acids and has a predicted molecular mass of 23.8 kDa. The apparent molecular mass of rh SLAMF7 is approximately 35-45 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human SLAMF7 / CRACC / CD319 Protein (His Tag) SDS-PAGE