

# Human Angiotensinogen/SerpinA8 (His Tag) recombinant protein



Catalog Number: 501652

## General Information

### Gene Name Synonym

Angiotensinogen (Serpine peptidase inhibitor, clade A, member 8), isoform CRA\_a

### Protein Construction

A DNA sequence encoding the human angiotensinogen (NP\_000020.1) (Met 1-Ala 485) was fused with a polyhistidine tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 98 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Asp 34

### Molecular Mass

The recombinant human angiotensinogen

comprises 463 amino acids and has a calculated molecular mass of 51.2 kDa. As a result of glycosylation, the apparent molecular mass of rh angiotensinogen is approximately 55-60 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human Angiotensinogen / SerpinA8 / AGT Protein (His Tag) SDS-PAGE