

H5N1 HA (His Tag) recombinant protein



Catalog Number: 501817

General Information

Protein Construction

A DNA sequence encoding the N-terminal segment (Met 1-Glu 339) of the influenza A hemagglutinin (A/Egypt/N05056/2009 (H5N1)) (ACT15357.1), termed as HA1, was fused with a C-terminal polyhistidine tag.

Organism

H5N1

Expression Host

Human Cells

QC Testing

Purity

> 98 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 17

Molecular Mass

The secreted recombinant HA1 subunit influenza A H5N1 (A/Egypt/N05056/2009 (H5N1)) HA comprises 334 amino acids and has a predicted

molecular mass of 38 kDa. As a result of glycosylation, it migrates as an approximately 45-50 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H5N1 (A/Egypt/N05056/2009) Hemagglutinin Protein (HA1 Subunit) (His Tag) SDS-PAGE