

Cynomolgus Fractalkine/CX3CL1 (His Tag) recombinant protein



Catalog Number: 501892

General Information

Protein Construction

A DNA sequence encoding the cynomolgus CX3CL1 (EHH60415.1) (Met1-Gly100) was expressed with a polyhistidine tag at the C-terminus. Cynomolgus and Rhesus CX3CL1 sequences are identical.

Organism

Cynomolgus

Expression Host

Human Cells

QC Testing

Purity

> (71.4+6.5+20.3) % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 25

Molecular Mass

The recombinant cynomolgus CX3CL1 consists of

87 amino acids and predicts a molecular mass of 10.1 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Cynomolgus / Rhesus Fractalkine / CX3CL1 Protein (His Tag)