

# Human IZUMO1 (His Tag) recombinant protein



Catalog Number: 501905

## General Information

### Gene Name Synonym

Oocyte binding/fusion factor; Sperm-specific protein izumo

### Protein Construction

A DNA sequence encoding the human IZUMO1 (AAH34769.1) (Met1-Arg292) was expressed with a polyhistidine tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Cys 22

### Molecular Mass

The recombinant human IZUMO1 consists of 282 amino acids and predicts a molecular mass of 32.1 KDa. It migrates as an approximately 43 KDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human IZUMO1 Protein (His Tag) SDS-PAGE