

General Information

Gene Name Synonym

GP protein; GP1,2; Virion spike glycoprotein

Protein Construction

A DNA sequence encoding the Zaire ebolavirus(strain Ebola virus H.sapiens-wt/SLE/2014/ManoRiver-G3686.1) GP (AIE11917.1) (Met1-Gln650) was expressed with a polyhistidine tag at the C-terminus.

Organism

EBOV

Expression Host

Human Cells

QC Testing

Purity

> (55.1+41.6) % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ile 33

Molecular Mass

The recombinant Zaire ebolavirus(strain Ebola virus H.sapiens-wt/SLE/2014/ManoRiver-G3686.1) GP consists 629 amino acids and predicts a molecular mass of 69.3 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Ebola virus EBOV (subtype Zaire, strain Ebola virus H.sapiens-wt/SLE/2014/ManoRiver-G3686.1) Glycoprotein / EBOV-G Protein (His Tag)