

Human ENO3 / beta-enolase (His Tag) recombinant protein



Catalog Number: 502206

General Information

Gene Name Synonym

2-phospho-D-glycerate hydro-lyase; Enolase 3;
Muscle-specific enolase; Skeletal muscle enolase

Protein Construction

A DNA sequence encoding the mature form of human ENO3 (AAH17249.1) (Met1-Lys434) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human ENO3 consists of 449 amino acids and predicts a molecular mass of 48.8 KDa. It migrates as an approximately 45 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 50mM Tris, 0.1% Brij35, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human ENO3 / beta-enolase Protein (His Tag)
SDS-PAGE