

Rat Neurexophilin-1/NXPH1 (His Tag) recombinant protein



Catalog Number: 502223

General Information

Protein Construction

A DNA sequence encoding the rat NXPH1 (Q63366-1) (Met 1-Gly 271) was expressed, fused with a polyhistidine tag at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 92 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 22

Molecular Mass

The secreted recombinant rat NXPH1 comprises 261 amino acids and predicts a molecular mass of 30 kDa. The apparent molecular mass of the rat

NXPH1 is approximately 45-50 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat Neurexophilin-1 / NXPH1 Protein (His Tag)
SDS-PAGE