

Human LRRN3 (His Tag) recombinant protein



Catalog Number: 502292

General Information

Gene Name Synonym

Neuronal leucine-rich repeat protein 3

Protein Construction

A DNA sequence encoding the human LRRN3 (AAH35133.1) extracellular domain (Met 1-Thr 628) was fused with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met 1

Molecular Mass

The recombinant human LRRN3 consists of 616 amino acids and predicts a molecular mass of 70 kDa as estimated in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 7.0, 10% gly

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human LRRN3 Protein (His Tag) SDS-PAGE