

General Information

Protein Construction

A DNA sequence encoding the human herpesvirus 4/Epstein-Barr virus gp350 (strain B95-8) (YP_401667.1) (Met1-Pro490) was expressed with a polyhistidine tag at the C-terminus.

Organism

EBV

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human herpesvirus 4/Epstein-Barr virus gp350 consists 501 amino acids and predicts a molecular mass of 53.7 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 % glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Epstein-Barr virus (Herpesvirus 4) EBV
Glycoprotein gp350 / EBV GP350 Protein (His Tag)