

General Information

Gene Name Synonym

Serpin A10

Protein Construction

A DNA sequence encoding the mature form of mouse SERPINA10 isoform 1 (Q8R121-1) (Met 1-Leu 448) was fused with a polyhistidine tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Phe 22

Molecular Mass

The secreted recombinant mouse SERPINA10

comprises 438 amino acids with a predicted molecular mass of 51 kDa. As a result of glycosylation, the apparent molecular mass of the protein is approximately 60-80 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse SerpinA10 / ZPI Protein (His Tag) SDS-PAGE