

Human TREML2/TLT-2 (His Tag) recombinant protein



Catalog Number: 502327

General Information

Gene Name Synonym

Triggering receptor expressed on myeloid cells-like protein 2

Protein Construction

A DNA sequence encoding the human TREML2 (NP_079083.2) extracellular domain (Met 1-Ser 268) was fused with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gly 19

Molecular Mass

The recombinant human TREML2 consists of 261

amino acids and has a predicted molecular mass of 28.5 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rh TREML2 is approximately 55-60 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human TREML2 / TLT2 Protein (His Tag) SDS-PAGE