

# Human BLMH/Bleomycin Hydrolase (His Tag) recombinant protein



Catalog Number: 502352

## General Information

### Protein Construction

A DNA sequence encoding the human BLMH (NP\_000377.1) (Ser 2-Glu 455) was expressed, with a polyhistide tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Activity

Measured by its ability to hydrolyze Met-AMC.  
The specific activity is >500 pmoles/min/μg.

### Purity

> 88 % as determined by SDS-PAGE

### Endotoxin

Please contact us for more information.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Met

### Molecular Mass

The recombinant human BLMH consisting of 461

amino acids and has a calculated molecular mass of 53.4 kDa as estimated in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human BLMH / BLM hydrolase Protein (His Tag)  
SDS-PAGE