

Human RUVBL1 / RVB1 (His Tag) recombinant protein



Catalog Number: 502397

General Information

Gene Name Synonym

49 kDa TATA box-binding protein-interacting protein; 54 kDa erythrocyte cytosolic protein; INO80 complex subunit H; Nuclear matrix protein 238; Pontin 52; TIP49a; TIP60-associated protein 54-alpha

Protein Construction

A DNA sequence encoding the human RUVBL1 (Q9Y265-1) (Met1-Lys456) was fused with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human RUVBL1 consists of 474 amino acids and has a calculated molecular mass of 52.4 kDa. The recombinant protein migrates as an approximately 57 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 7.4, 10% glycerol

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human RUVBL1 / RVB1 Protein (His Tag) SDS-PAGE