

H7N7 HA (His Tag) recombinant protein



Catalog Number: 502463

General Information

Protein Construction

A DNA sequence encoding the N-terminal segment (Met 8-Arg 345) of Influenza A virus (A/Netherlands/219/03 (H7N7)) hemagglutinin (AAR02640.1), termed as HA1, was fused with a polyhistidine tag at the C-terminus.

Organism

H7N7

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 26

Molecular Mass

The recombinant HA1 subunit of Influenza A virus hemagglutinin (A/Netherlands/219/03 (H7N7))

comprises 331 amino acids with the predicted molecular mass of 36.4 kDa. As a result of glycosylation, it migrates as an approximately 48 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H7N7 (A/Netherlands/219/03) Hemagglutinin / HA Protein (His Tag) SDS-PAGE