

# Mouse CD55/DAF (His Tag) recombinant protein



Catalog Number: 502468

## General Information

### Protein Construction

A DNA sequence encoding the mature form of mouse CD55 (Q61475) (Met 1-Thr 361) was expressed, with a C-terminal polyhistidine tag.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

> 92 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Asp 35

### Molecular Mass

The secreted recombinant mouse CD55 comprises 338 amino acids and has a calculated molecular mass of 37.2 kDa. As a result of glycosylation, the

apparent molecular mass of rmCD55 is approximately 60 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse CD55 / DAF Protein (His Tag) SDS-PAGE