

Human IFNA10 (Fc Tag) recombinant protein



Catalog Number: 502587

General Information

Gene Name Synonym

Interferon alpha-6L; Interferon alpha-C

Protein Construction

A DNA sequence encoding the human IFNA10 (P01566) (Cys 24-Asp 189) was fused with the Fc region of human IgG1 at the N-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Activity

Measured in antiviral assays using WISH human amnion cells infected with vesicular stomatitisvirus (VSV).

The ED₅₀ for this effect is 20-80 pg/mL.

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu

Molecular Mass

The recombinant human IFNA10/Fc chimera is a disulfide-linked homodimer. The reduced monomer comprises 427 amino acids with a predicted molecular mass of 48 kDa.

AsrhIFNA10/Fc monomer migrates as an approximately 50 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C.

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human Interferon alpha 10 / IFNA10 Protein (Fc Tag) SDS-PAGE