

Mouse GAPDH (His Tag) recombinant protein



Catalog Number: 502658

General Information

Gene Name Synonym

Peptidyl-cysteine S-nitrosylase GAPDH

Protein Construction

A DNA sequence encoding the mouse GAPDH (P16858) (Met1-Glu333) was expressed with a polyhistidine tag at the N-terminus.

Organism

Mouse

Expression Host

E. coli

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant mouse GAPDH consists of 349 amino acids and predicts a molecular mass of 37.9

KDa. It migrates as an approximately 38 KDa band in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile 50 mM Tris, 200 mM Arg, 10% gLycerol.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse GAPDH Protein (His Tag) SDS-PAGE