

Human Discoidin / DCBLD1 (His Tag) recombinant protein



Catalog Number: 502746

General Information

Protein Construction

A DNA sequence encoding the human DCBLD1 (Q8N8Z6-2) (Met1-Val459) was expressed with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu 35

Molecular Mass

The recombinant human DCBLD1 consists of 436 amino acids and predicts a molecular mass of 48.4 KDa. It migrates as an approximately 65 KDa band

in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human Discoidin / DCBLD1 Protein (His Tag)
SDS-PAGE