

Human DBH / Dopamine beta-Hydroxylase (His Tag) recombinant protein



Catalog Number: 502758

General Information

Gene Name Synonym

Dopamine beta-monooxygenase; Soluble dopamine beta-hydroxylase

Protein Construction

A DNA sequence encoding the human DBH (P09172) (Ser26-Gly603) was expressed with an N-terminal polyhistidine tag.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human DBH comprises 598 amino acids and has a predicted molecular mass of 67.3 kDa. The apparent molecular mass of the protein is approximately 68 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human DBH / Dopamine beta-Hydroxylase Protein (His Tag) SDS-PAGE