

Mouse CADM4/IGSF4C/NECL-4 (His Tag) recombinant protein



Catalog Number: 502842

General Information

Gene Name Synonym

Immunoglobulin superfamily member 4C; Nectin-like protein 4; TSLC1-like protein 2

Protein Construction

A DNA sequence encoding the extracellular domain of mouse CADM4 (NP_694752.1) (Met 1- Ala 324) was expressed, with a polyhistidine tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu 26

Molecular Mass

The recombinant mouse CADM4 consists of 310

amino acids and has a predicted molecular mass of 34.5 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rm CADM4 is approximately 50-60 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse CADM4 / IGSF4C / TSL2 Protein (His Tag) SDS-PAGE