

# Human UBE2A/HHR6A (His Tag) recombinant protein



Catalog Number: 502993

## General Information

### Gene Name Synonym

RAD6 homolog A; Ubiquitin carrier protein A; Ubiquitin-protein ligase A

### Protein Construction

A DNA sequence encoding the mature form of human UBE2A (P49459) (Met 1-Cys 152) was expressed, with a polyhistidine tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

> 80 % as determined by SDS-PAGE

### Endotoxin

Please contact us for more information.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Met 1

### Molecular Mass

The recombinant human UBE2A consisting of 167 amino acids and has a calculated molecular mass of 19.2 kDa. It migrates as an approximately 18.5 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, 20% glycerol, pH 7.5  
1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human UBE2A / HHR6A Protein (His Tag) SDS-PAGE