

Mouse FAM3D / Oit1 (His Tag) recombinant protein



Catalog Number: 503000

General Information

Gene Name Synonym

Oncoprotein-induced protein 1; Protein EF-7

Protein Construction

A DNA sequence encoding the mouse FAM3D (P97805) (Met1-Met223) was expressed with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Thr 26

Molecular Mass

The recombinant mouse FAM3D comprises 209 amino acids and has a predicted molecular mass of 23.5 kDa. The apparent molecular mass of the protein is approximately 23 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse FAM3D / Oit1 Protein (His Tag) SDS-PAGE