

H2N2 HA (His Tag) recombinant protein



Catalog Number: 503059

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/Guiyang/1/1957(H2N2)) hemagglutinin (ACD85231.1) (Met1-Arg340), termed as HA1, was expressed with a C-terminal polyhistidine tag.

Organism

H2N2

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 16

Molecular Mass

The recombinant HA1 subunit of Influenza A virus (A/Guiyang/1/1957(H2N2)) comprises 336 amino acids and has a predicted molecular mass of 37.8 kDa. The apparent molecular mass of the protein

is approximately 49.2 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H2N2 (A/Guiyang/1/1957)
Hemagglutinin Protein (HA1 Subunit) (His Tag)
SDS-PAGE