

## General Information

### Gene Name Synonym

Fatty acid-binding protein 6; Ileal lipid-binding protein; Intestinal 15 kDa protein; Intestinal bile acid-binding protein

### Protein Construction

A DNA sequence encoding the human FABP6 (AAH22489.1) (Met 1-Ala 128) was expressed and purified.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

> 98 % as determined by SDS-PAGE

### Endotoxin

Please contact us for more information.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Met 1

### Molecular Mass

The recombinant human FABP6 consisting of 128 amino acids and has a calculated molecular mass of 14.4 kDa as estimated in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.5

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human FABP6 / I-BABP Protein SDS-PAGE