

General Information

Gene Name Synonym

GP1,2; GP1; GP2; GP2-delta

Protein Construction

A DNA sequence encoding the Zaire ebolavirus (strain Mayinga 1976) GP (AAC54887.1) (Met1-Phe308) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

EBOV

Expression Host

Human Cells

QC Testing

Purity

> (77.5+ 19.9) % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ile 33

Molecular Mass

The recombinant Zaire ebolavirus (strain Mayinga 1976) GP consists 514 amino acids and predicts a molecular mass of 57.6 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.