

# H13N8 HA (His Tag) recombinant protein



Catalog Number: 503152

## General Information

### Protein Construction

A DNA sequence encoding the Influenza A virus (A/black-headed gull/Netherlands/1/00 (H13N8)) hemagglutinin (AAV91212.1) (Met1-Gln528), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

### Organism

H13N8

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Activity

Measured by its ability to agglutinate guinea pig red blood cells.

HA titer is 0.05-0.25 $\mu$ g/mL for 1%GRBC

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Asp 19

### Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/black-headed gull/Netherlands/1/00 (H13N8)) consists 521 amino acids and predicts a molecular mass of 58.5 kDa.

### Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 % glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Influenza A H13N8 (A/black-headed gull/Netherlands/1/00) Hemagglutinin / HA Protein (His Tag)