

H4N8 HA (His Tag) recombinant protein



Catalog Number: 503225

General Information

Gene Name Synonym

Hemagglutinin HA1 chain; Hemagglutinin HA2 chain

Protein Construction

A DNA sequence encoding the extracellular domain of influenza A virus hemagglutinin (A/chicken/Alabama/1/1975 (H4N8)) (P19695) (Met 1-Ile528) was fused with a polyhistidine tag at the C-terminus.

Organism

H4N8

Expression Host

Human Cells

QC Testing

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 17

Molecular Mass

The secreted recombinant hemagglutinin of

Influenza A virus H4N8 HA

(A/chicken/Alabama/1/1975(H4N8)) comprises 523 amino acids and has a predicted molecular mass of 59 kDa. As a result of glycosylation, it migrates as an approximately 75 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H4N8 (A/chicken/Alabama/1/1975) Hemagglutinin / HA Protein (His Tag) SDS-PAGE