

## General Information

### Gene Name Synonym

Interleukin 2 receptor, gamma (Severe combined immunodeficiency), isoform CRA\_a; Protein IL2rg

### Protein Construction

A DNA sequence encoding the rat IL2RG (Q68FU6) (Met1-Ala262) was expressed, fused with a polyhistidine tag at the C-terminus.

### Organism

Rat

### Expression Host

Human Cells

## QC Testing

### Purity

> 99 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Ser 25

### Molecular Mass

The recombinant rat IL2RG comprises 249 amino acids and predicts a molecular mass of 29.5 kDa. The apparent molecular mass of the recombinant protein is approximately 41-50 kDa in SDS-PAGE under reducing conditions due to glycosylation.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Rat IL2RG / CD132 Protein (His Tag) SDS-PAGE