

Canine IGFBP5 / IGFBP-5 (Fc Tag) recombinant protein



Catalog Number: 503256

General Information

Protein Construction

A DNA sequence encoding the canine IGFBP5 (F1PEQ5) (Met1-Glu271) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Canine

Expression Host

Human Cells

QC Testing

Purity

(26.4+14.4+6.3+4.6+46.8) % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Leu 20

Molecular Mass

The recombinant canine IGFBP5/Fc is a disulfide-linked homodimer. The reduced monomer comprises 493 amino acids and has a predicted

molecular mass of 55.5 kDa. The apparent molecular mass of the protein is approximately 64, 45, 43, 40, and 33 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Canine IGFBP5 / IGFBP-5 Protein (Fc Tag) SDS-PAGE