

# Mouse IGSF8 (His Tag) recombinant protein



Catalog Number: 503368

## General Information

### Gene Name Synonym

CD81 partner 3; Glu-Trp-Ile EWI motif-containing protein 2; Keratinocyte-associated transmembrane protein 4; Prostaglandin regulatory-like protein

### Protein Construction

A DNA sequence encoding the mouse IGSF8 (AAH48387.1) (Met1-Thr577) was expressed with a C-terminal polyhistidine tag.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Arg 26

### Molecular Mass

The recombinant mouse IGSF8 comprises 563 amino acids and has a predicted molecular mass of 60.2 kDa. The apparent molecular mass of the protein is approximately 66 kDa in SDS-PAGE under reducing conditions due to glycosylation.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse PGRL / IGSF8 Protein (His Tag) SDS-PAGE