

Human CD32b/Fc gamma RIIB/FCGR2B (His & AVI Tag) recombinant protein



Catalog Number: 503484

General Information

Gene Name Synonym

CDw32; Fc-gamma RII-b

Protein Construction

A DNA sequence encoding the extracellular domain (Ala 46-Ile 224) of mature human CD32b (NP_001002274.1) protein was fused with a signal peptide at the N-terminus and a c-terminal polyhistidine tagged AVI tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Activity

Immobilized human CD32b-AVI-His (Cat:503484) at 10 µg/ml (100 µl/well) can bind biotinylated Human IgG1, The EC₅₀ of biotinylated Human IgG1 is 0.20-0.48 µg/ml.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 46

Molecular Mass

The recombinant mature human CD32b consists of 214 amino acids after proteolytic removal of the signal peptide and predicts a molecular mass of 24 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of the protein is approximately 30-35 kDa protein due to glycosylation.

Formulation

- Lyophilized from sterile PBS, pH 7.4
1. 5 % trehalose and mannitol are added as protectants before lyophilization.
 2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human CD32b / FCGR2B Protein (His & AVI Tag)
SDS-PAGE