

General Information

Gene Name Synonym

CD218 antigen-like family member A; CDw218a; IL1 receptor-related protein

Protein Construction

A DNA sequence encoding the extracellular domain of human IL18R1 (NP_003846.1) (Met 1-Arg 329) was expressed, fused with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 19

Molecular Mass

The recombinant human IL18R1 comprises 322

amino acids with a predicted molecular mass of 37.2 kDa. As a result of glycosylation, the apparent molecular mass of rhIL18R1 is approximately 60-70 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human IL18R1 / CD218a Protein (His Tag) SDS-PAGE