

H1N1 HA (His Tag) recombinant protein



Catalog Number: 503682

General Information

Gene Name Synonym

Hemagglutinin HA1 chain; Hemagglutinin HA2 chain

Protein Construction

A DNA sequence encoding an N-terminal segment of Influenza A virus H1N1 hemagglutinin (A/Brevig Mission/1/1918 (H1N1)) (AAD17229.1) (Met 1-Arg 344), termed as HA1 subunit, was fused with a polyhistidine tag at the C-terminus.

Organism

H1N1

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 18

Molecular Mass

Recombinant HA1 subunit of Influenza A virus

H1N1 hemagglutinin (A/Brevig Mission/1/1918 (H1N1)) comprises 338 amino acids with the predicted molecular mass of 37.5 kDa. As a result of glycosylation, it migrates as an approximately 53 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4
1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H1N1 (A/Brevig Mission/1/1918) Hemagglutinin Protein (HA1 Subunit) (His Tag) SDS-PAGE