

Cynomolgus FGFR3/CD333 (Fc Tag) recombinant protein



Catalog Number: 503694

General Information

Protein Construction

A DNA sequence encoding the cynomolgus FGFR3 (XP_001101108.1) (Met1-Gly375) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Cynomolgus

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu 23

Molecular Mass

The recombinant cynomolgus FGFR3 is a disulfide-linked homodimer. The reduced monomer comprises 594 amino acids and has a

calculated molecular mass of 65.1 KDa. The apparent molecular mass of the protein is approximately 99.9 KDa respectively in SDS-PAGE.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Cynomolgus FGFR3 Protein (Fc Tag) SDS-PAGE