

## General Information

### Protein Construction

A DNA sequence encoding the HA (CAC84244.1) (Met1-Gln529) was expressed with a polyhistidine tag at the C-terminus.

### Organism

H6N4

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Activity

Measured by its ability to agglutinate guinea pig red blood cells.

HA titer is 0.5-2 $\mu$ g/mL for 1%GRBC

### Purity

> 90 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Ile 17

### Molecular Mass

The recombinant HA consists 464 amino acids and predicts a molecular mass of 52.6 kDa.

### Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 %glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Influenza A H6N4 (A/chicken/Hong Kong/17/1977)  
Hemagglutinin / HA Protein (His Tag)