

# H3N2 Neuraminidase (Active) recombinant protein



Catalog Number: 503738

---

## General Information

### Protein Construction

A DNA sequence encoding the Influenza A virus (A/Babol/36/2005 (H3N2)) neuraminidase (ACN50232.1) (His 36-Pro 459) was expressed, the cell lysates are collected, and bio-activity was tested. There is an amino acid change from Asparagine to Serine (N294S mutation) in NA / Neuraminidase.

### Organism

H3N2

### Expression Host

Human Cells

## QC Testing

### Activity

Measured by its ability to cleave a fluorogenic substrate, 2'-(4-Methylumbelliferyl)- $\alpha$ -D-N-acetylneuraminic acid.

The specific activity is > 80 U

The specific activity is > 500 U

The specific activity is > 800 U

One unit is defined as the amount of enzyme required to cleave 1 nmole of 2'-(4-Methylumbelliferyl)- $\alpha$ -D-N-acetylneuraminic acid per minute at pH 7.5 at 37°C.

### Endotoxin

< 1.0 EU per  $\mu$ g of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Molecular Mass

The influenza H3N2 virus Neuraminidase comprises 443 amino acids.

### Formulation

Lyophilized from sterile PBS, 0.6% Triton X-100, 7% Trehalose, 6% Mannitol, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

It is recommended that 1 ml sterile water be added to the vial to prepare a stock solution.