

Human EphB4/Eph Receptor B4 (aa 563-987, His & GST Tag) recombinant protein



Catalog Number: 503772

General Information

Gene Name Synonym

Hepatoma transmembrane kinase; Tyrosine-protein kinase TYRO11

Protein Construction

A DNA sequence encoding the human EPHB4 (P54760) (Leu563-Tyr987) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

The specific activity was determined to be 47 nmol/min/mg using Poly(Glu:Tyr) 4:1 as substrate.

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human EPHB4/GST chimera consists of 662 amino acids and has a calculated molecular mass of 75.2 kDa. The recombinant protein migrates as an approximately 66 kDa band in SDS-PAGE under reducing conditions.

Formulation

Supplied as sterile 20mM Tris, 500mM NaCl, pH 8.0, 3mM DTT, 10% gly

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human EphB4 / HTK Protein (aa 563-987, His & GST Tag) SDS-PAGE