

General Information

Protein Construction

A DNA sequence encoding the Influenza B virus (B/Massachusetts/03/2010) hemagglutinin (AGB07680.1) (Met1-Thr548), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

Organism

Influenza B

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

Measured by its ability to agglutinate guinea pig red blood cells.

HA titer is 0.01-0.05 μ g/mL for 1%GRBC

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μ g protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 16

Molecular Mass

The recombinant hemagglutinin of Influenza B virus (B/Massachusetts/03/2010) consists 544 amino acids and predicts a molecular mass of 59.1 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, pH 7.0, 10 % glycerol.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza B (B/Massachusetts/03/2010)
Hemagglutinin / HA Protein (His Tag)