

H10N3 HA (His Tag) recombinant protein



Catalog Number: 503818

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/duck/Hong Kong/786/1979(H10N3)) hemagglutinin (BAF46762.1) (Met1-Ile525) termed as HA, was expressed with a polyhistidine tag at the C-terminus.

Organism

H10N3

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Leu 17

Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/duck/Hong Kong/786/1979(H10N3))

consists 520 amino acids and predicts a molecular mass of 57.9 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 % glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H10N3 (A/duck/Hong Kong/786/1979) Hemagglutinin / HA Protein (His Tag)