

Mouse LAIR-1 / CD305 (Fc Tag) recombinant protein



Catalog Number: 503921

General Information

Protein Construction

A DNA sequence encoding the mouse LAIR1 (Q8BG84-6) (Gln22-Ser133) was expressed, fused with the Fc region of human IgG1 at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 22

Molecular Mass

The recombinant mouse LAIR1 /Fc is a disulfide-linked homodimer. The reduced monomer comprises 353 amino acids and has a predicted

molecular mass of 39.6 kDa. The apparent molecular mass of the protein is approximately 49-53 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse LAIR1 Protein (Fc Tag) SDS-PAGE