

Canine CD40/TNFRSF5 (His Tag) recombinant protein



Catalog Number: 503987

General Information

Gene Name Synonym

B-cell surface antigen CD40; CD40L receptor

Protein Construction

A DNA sequence encoding the canine CD40 (Q7YRL5) (Met1-Ala194) was expressed with a C-terminal polyhistidine tag.

Organism

Canine

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu 21

Molecular Mass

The recombinant canine CD40 comprises 185

amino acids and has a predicted molecular mass of 20.5 kDa. The apparent molecular mass of the protein is approximately 33 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Canine CD40 / TNFRSF5 Protein (His Tag) SDS-PAGE