

Human EphB4/Eph Receptor B4 recombinant protein



Catalog Number: 504023

General Information

Gene Name Synonym

Hepatoma transmembrane kinase; Tyrosine-protein kinase TYRO11

Protein Construction

The mature form of human EphB4 (NP_004435.3) extracellular domain (Met 1-Ala 539) with quinary-aa peptide (DDDDK) at the C-terminus was expressed and purified.

Organism

Human

Expression Host

Human Cells

QC Testing

Activity

Measured by its binding ability in a functional ELISA. Immobilized human EphB4 at 2 µg/ml (100 µl/well) can bind human EphrinB2 with a linear range of 1-125 ng/ml.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Leu 16

Molecular Mass

The recombinant human EphB4 is a monomeric protein consisting of 530 amino acids and has a calculated molecular mass of 57.7 kDa. As a result of glycosylation, rhEphB4 migrates as an approximately 65-70 kDa protein in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 100mM NaCl, 50mM Tris, pH 7.5

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human EphB4 / HTK Protein SDS-PAGE