

# Canine IL13RA2 / CD213A2 (His Tag) recombinant protein



Catalog Number: 504083

## General Information

### Protein Construction

A DNA sequence encoding the canine IL13RA2 (Q95LF0) (Met1-Thr338) was expressed with a C-terminal polyhistidine tag.

### Organism

Canine

### Expression Host

Human Cells

## QC Testing

### Activity

Measured by its ability to inhibit IL13 dependent proliferation of TF-1 human erythroleukemic cells. The ED<sub>50</sub> for this effect is typically 0.01-0.04 µg/mL.

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Ser 22

### Molecular Mass

The recombinant canine IL13RA2 comprises 328 amino acids and has a predicted molecular mass of 38.6 kDa. The apparent molecular mass of the protein is approximately 39-46 kDa in SDS-PAGE under reducing conditions due to glycosylation.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Canine IL13RA2 / IL13R Protein (His Tag) SDS-PAGE