

Rat CDH13 / Cadherin-13 / H Cadherin recombinant protein



Catalog Number: 504213

General Information

Gene Name Synonym

Protein Cdh13; T-cadherin

Protein Construction

A DNA sequence encoding the rat CDH13 (NP_620244.1)(Met1-Ala692) was expressed with six amino acids (LEVLFQ) at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 22

Molecular Mass

The recombinant rat CDH13 comprises 678 amino

acids and predicts a molecular mass of 74.3 kDa. The apparent molecular mass of the recombinant protein is approximately 84-96 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat CDH13 / Cadherin-13 / H Cadherin Protein
SDS-PAGE