

H1N1 HA (His Tag) recombinant protein



Catalog Number: 504235

General Information

Gene Name Synonym

Hemagglutinin HA1 chain; Hemagglutinin HA2 chain

Protein Construction

A DNA sequence encoding the Influenza A virus (A/Puerto Rico/8/1934(H1N1)) hemagglutinin (P03452.2) (Met1-Gln528), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

Organism

H1N1

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

Measured by its ability to agglutinate guinea pig red blood cells.

HA titer is 0.2-1 μ g/mL for 1%GRBC

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μ g protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 18

Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/Puerto Rico/8/1934(H1N1)) consists 522 amino acids and predicts a molecular mass of 59.2 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 % glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H1N1 (A/Puerto Rico/8/1934)
Hemagglutinin / HA Protein (His Tag)